

The Radixact[®] System

High-Fidelity images
Helical fan-beam MVCT
and kVCT

Longest treatment fields
Continuous couch movement
with fully rotating gantry

**Superior control of
dose delivery**
Dynamic jaws and fast
binary MLC

Quick installation
RF chain pre-calibrated
and ships intact



Fast and flexible imaging
Collimated kV source and
movable flat panel detector

Fast treatment delivery
1,000 MU/minute linac

**Dose minimization to
organs at risk**
True helical delivery from
any angle with continuous
couch movement

Tight personalized margins
ClearRT[™] kVCT imaging and
Synchrony[®] real-time
delivery adaptation

True Helical-Based,
High-Fidelity
Imaging and Ultra-
Precise Delivery